

**IgG Fraction of Rabbit Anti-M<sub>1</sub> Receptor Serum**

AS-3701G

Lot # 8360

This IgG fraction was prepared by precipitation, dialysis and chromatography from an antiserum that was raised in a rabbit which was immunized with a peptide analogue of the carboxyl terminal of the M<sub>1</sub> receptor covalently attached onto a carrier protein. The antiserum has been shown to be specific for the COOH terminal of the M<sub>1</sub> receptor and has been characterized by western blotting, ELISA, and cell staining techniques. This IgG fraction is suitable for immunocytochemical and western immunoblotting detection of the receptor. Dilute the lyophilized antibody with 0.1 ml of 10 mg/ml BSA in PBS for the equivalent concentration in whole antiserum, or with additional buffer for more dilute antibody. This vial contains 1.58 mg of purified IgG.

**Antiserum Specificity**

<b>Polypeptide</b>	<b>% Cross Reactivity</b>
M1 Muscarinic Receptor (451-460)	100
M1 Muscarinic Receptor	60
M2 Muscarinic Receptor (457-466)	0
M2 Muscarinic Receptor	0
M3 Muscarinic Receptor (580-589)	0
M3 Muscarinic Receptor	0
M4 Muscarinic Receptor (469-478)	0
M4 Muscarinic Receptor	0
M5 Muscarinic Receptor (519-531)	0
M5 Muscarinic Receptor	0

**Immunocytochemical Staining**

This IgG fraction has been found to stain PLP fixed rat brain sections known to express the muscarinic M<sub>1</sub> receptor using ABC techniques at a concentration of 5 µg/ml. For information on PLP fixative see the FAQ page on our web site at [www.RDAbs.com](http://www.RDAbs.com).

**Western Immunoblot**

Western immunoblots using whole rat brain homogenate have been successful at a concentration of 20 µg/ml.

**Western Blotting Protocol**

1. After SDS-PAGE (on either 4-15% gradient gels or single percentage gels, such as 12% gels) and electrophoretic transfer to PVDF, block the membrane overnight with 4% normal goat serum using TBS/Tween-20 buffer as diluent.
2. Wash x 2 with TBS/Tween-20.
3. Apply the IgG fraction after diluting to at least 20 µg/ml (Note: higher dilutions may be needed). Use 2% normal goat serum in TBS/Tween-20 as buffer for the primary antibody. Let the primary antibody bind for 2-4 hours.
4. Wash x 3 with TBS/Tween-20.
5. Apply affinity purified HRP-goat anti-rabbit IgG antiserum diluted 1:2500 (dilution may vary depending upon supplier) in 2% normal goat serum in TBS/Tween-20. Incubate 1-2 hours. Note: greater sensitivity may be achieved using ABC techniques.
6. Wash x 4 for 5 minute per wash cycle in TBS/Tween-20.
7. Develop color using the enhanced DAB reaction.