

Rabbit Anti-Calcitonin (Salmon) Sera

AS-1218

Lot # 4597

These antisera were raised in rabbits which were immunized with Salmon Calcitonin covalently attached onto a carrier protein. The antisera have been shown to be specific for the immunocytochemical detection of Calcitonin. This antisera has been characterized by ELISA techniques. Dilute the lyophilized antibody with 0.1 ml of 10 mg/ml BSA in PBS for the equivalent of whole antiserum, or with additional buffer for more dilute antibody.

Antiserum Specificity

Polypeptide	% Cross Reactivity
Calcitonin (Salmon)	100
Calcitonin (Human)	0
CGRP (Rat)	0
ANP	0
PTH	0
Somatostatin 28	0
Vasoactive Intestinal Peptide	0